


Transgene	CAG promoter
-----------	--------------

primer	Sequence(5'-3')	Size	target gene	accession
CAGGS-F CAGGS-R	TAATCAATTACGGGGTCATTAGTTCATAGC TCCATAAGGTCATGTACTGGGCATAATGC	301bp	pCAGGS	
comment				

PCR condition						
Taq polymerase	BIOTAQ™ HS DNA Polymerase(BIOLINE, London, UK)					
Thermal cycler	PC-808(ASTEC)					
PCR buffer	Ampdirect Plus					
	first denature	denature	anneal	extension	cycle	final extension
PCR	°C	94	94	55	35	72
	min	3	0.5	1		0.75
comment	PCR method; Blood was applied to FTA®card (GE Healthcare UK Ltd.,UK) and dried. 1.5mm FTA disc was removed from the bloody-stained region on FTA®card. Untreated sample discs were placed directly in 15 μ L PCR mixture containig 1 × Ampdirect®Plus, 0.2 μ M each primer and 0.4 units of BIOTAQ™ HS DNA Polymerase.					

Electrocataphoresis ;4% agarose gel 300V, 30min	
	
Strain	LEW-Tg(CAG-EGFP)1Ys W-Tg(CAG-cre)81Jmsk