

| | | | | |
|--------------------------|--|-------|-------------|-------------|
| Gene symbol | <i>Atrn</i> ^{mv} | | | |
| primer | Sequence(5'-3') | Size | target gene | accession |
| EST13-117 rAtrn-mv-12 | TGCATCCTGTAAGCCAATT TGCAGGTGCAGGTAAGACAG | 528bp | <i>Atrn</i> | NM_031351.1 |
| EST13-117 EST13-113 | TGCATCCTGTAAGCCAATT GCAGGCCAGGTCTATGTGTA | 576bp | <i>Atrn</i> | NM_031351.1 |
| comment | EST13-117&rAtrn-mv-26はミュータントアレルの検出 EST13-117&EST13-113は野生型アレルの検出 exon1を含む6914bpの欠失 | | | |

| PCR condition | | | | | | |
|----------------|-----------------------------------|----------------|----------|--------|-----------|-------|
| Taq polymerase | Nova Tag Hot start DNA polymerase | | | | | |
| Thermal cycler | PC-808(ASTEC) | | | | | |
| PCR buffer | Ampdirect Plus | | | | | |
| | | first denature | denature | anneal | extension | cycle |
| PCR | °C | 94 | 94 | 60 | 72 | 35 |
| | min | 3 | 0.5 | 1 | 0.75 | |
| comment | FTAカードに採取した血液を未精製のままテンプレートとして用いる。 | | | | | |

| | | |
|---|-----------|-------------|
| Electrocataphoresis ;4% agarose gel 300V, 30min | | |
| | | |
| Primer | EST13-117 | EST13-117 |
| | EST13-113 | rAtrn-mv-12 |
| Size | 576 bp | 528 bp |
| mv/mv | + | - |
| mv/+ | + | + |
| +/+ | - | + |
| Strain harboring mutation | MV/Opu | |

